


# ***Water Distilling***

---

## **Model OPTI-M**



### **Glass Distillation with Metal Heater**

- ✓ Mounted on a sturdy metal stand and output 4 l/h
- ✓ Temperature : 30 - 40°C, which is taken care by efficient cooling by coiled Borosilicate condenser
- ✓ Conductivity of water: 1.5 to 2.5M, PH 5 - 6
- ✓ Pyrogen free, Power consumption 3kW Single Phase through a metal heater (chromium plated), which is inside a horizontal Glass Boiler.