



Vertical Laminar Flow Cabinet



IBS-V800



IBS-DDC / IBS-SDC

TECHNICAL PARAMETERS :

Model	IBS-V800	IBS-DDC	IBS-SDC
External Size (WxDxH), mm	802x650x1550	1040x615x1770	1440x615x1770
Internal Size (WxDxH), mm	800x530x540	940x540x545	1340x540x545
Work Surface Height	660mm	750mm	750mm
Display	LED Display		
Airflow Velocity	Average of 0.3~0.5m/s, air speed adjustable.		
Material	Main Body : Cold-rolled steel with anti-bacteria powder coating		
	Work Table : 304 stainless steel		
Pre-filter	Polyester fiber, washable		
HEPA Filter	99.999% efficiency at 0.3um		
Noise	<60dB		
Front Window	Manual,5mm toughened glass, anti-UV		
Max Opening	490mm	310mm	310mm
Illuminating Lamp	LED Lamp 8Wx1	LED Lamp 12Wx1	LED Lamp 16Wx1
	20Wx1	18Wx1	30Wx1
UV Lamp	Emission of 253.7 nanometers		
Consumption	350W	350W	600W
Standard Accessory	LED Lamp	LED Lamp	LED Lamp
	UV lampx2	UV lampx2	UV lampx2
	Base stand	Base stand	Base stand
Optional Accessory	Electric height adjustable bass stand		
Caster	Universal wheel with leveling feet		
Power Supply	AC 220V±10%, 50/60 Hz; 110V±10%, 60HZ		
Gross Weight	116kg	131kg	174kg
Package Size (WxDxH)mm	960x800x1800	1200x850x1360	1600x850x1360