



# Vertical Laminar Flow Cabinet-Double Sides Type

## ADVANTAGE

1. With memory function in case of power-failure.
2. Double sides work bench, 2~4 operators can work together.



**LED Display**  
Microprocessor control system, LED display



**UV Lamp**  
Emission of 253.7 nanometers for most efficient decontamination



**Double Sides**  
Double sides work bench, operators can sit face to face



## TECHNICAL PARAMETERS :

Model	IBS-DSC	BBS-SSC
External Size (WxDxH)	1040x660x1770mm	1440x660x1770mm
Internal Size (WxDxH)	940x560x545mm	1340x560x545mm
Work Surface Height	750mm	
Display	LED Display	
Airflow Velocity	Average of 0.3-0.5m/s	
Material	Main Body : Cold-rolled steel with anti-bacteria powder coating	
	Work Table: 304 stainless steel	
Pre-filter	Polyester fiber, washable	
HEPA Filter	99.999% efficiency at 0.3um	
Noise	<60dB	
Front Window	Manual, 5mm toughened glass, anti-UV	
Max Opening	310mm	310mm
UV Lamp	18Wx1	30Wx1
	Emission of 253.7 nanometers	
LED Lamp	12Wx1	16Wx1
Consumption	350W	600W
Caster	Universal caster with leveling feet	
Power Supply	AC 220V±10%, 50/60 Hz; 110V±10%, 60 HZ	
Standard Accessory	LED lamp, UV lampx2, Base stand	
Optional Accessory	Electric height adjustable base stand	
Gross Weight	150kg	186kg
Package Size (WxDxH)	1190x890x1340mm	1590x890x1340mm