



Laboratory Equipment

TDS Meter WF139A

- Measuring range : 100-19900 $\mu\text{S}/\text{cm}$
- Accuracy : $\pm 2\%$ F·S
- Resolution : 1ppm (mg/L)
- Operating Temperature : 0-50°C
- Power Supply : battery 4×1.5V (AG13)
- Dimensions : 165×26×26mm
- Weight : 63g



INNOTECH Inc

2648E Workman Ave # 413

West Covina, CA 91791

Tel. : 626-420759

Email : innotech@innotechlab.com

<http://www.innotechlab.com>