



Laboratory Equipment

TDS Meter WF139

- Measuring range : 100-19900 μ S/cm
- Accuracy : \pm 2% F·S
- Resolution : 1ppm (mg/L)
- Operating Temperature : 0-50 $^{\circ}$ C
- Temperature Compensation : 0-50 $^{\circ}$ C
- Power Supply 2 \times 3V (CR2032)
- Dimensions : 150 \times 27 \times 20mm
- Weight : 46g



INNOTECH Inc

2648E Workman Ave # 413

West Covina, CA 91791

Tel. : 626-420759

Email : innotech@innotechlab.com

<http://www.innotechlab.com>