

Heating and Drying Oven



Features:

1. Microcomputer control (PID controlling). LED display, more precision and stability. Audio and visual alarm. Two stainless steel shelves. Height adjustable.
2. External Material : excellent stainless steel panel, round angle structure, smooth and easy to clean.
3. Advanced air passage, low noise, stable performance, adjustable airflow.

Model	BJPX-Spring	BJPX-Summer	BJPX-Autumn	BJPX-Winter
Temperature range	Chamber temperature +5°C - 300°C			
Temperature control	0.1°C			
Temperature fluctuation	±0.5°C			
Shelves	2	2	2	2
Consumption power	1600W	2200W	2400W	4800W
Power	110V-220V / 50-60Hz			
Working time	Continuous			
Chamber size (mm)	400 x 400 x 450	450 x 550 x 550	500 x 600 x 750	800 x 700 x 1140
External dimension (mm)	590 x 540 x 790	640 x 690 x 890	690 x 740 x 1090	1082 x 977 x 1330



INNOTECH Inc.

2648E Workman Ave # 413

West Covina, CA 91791

Tel : 626-420756

Email : innotech@innotechlab.com