



Laboratory Equipment

pH Meter WF03II

- Measuring range: 0.00~14.00 pH
- Accuracy ± 0.1 pH
- Resolution 0.01 pH
- Operating Temperature 0~50°C
- Temperature Compensation 0~50°C
- Power Supply 4×1.5V (AG13)
- Dimensions: 165×29×26mm



INNOTECH Inc

2648E Workman Ave # 413

West Covina, CA 91791

Tel. : 626-420759

Email : innotech@innotechlab.com

<http://www.innotechlab.com>