



Refractometer



- **Model REF-101**
Portable Refractometer with 0 to 10% Brix measuring range
- **Model VBR32**
Handheld Refractometer with 0~32% Brix measuring range, 0.2% resolution
- **Model REF-104**
Portable Refractometer with 28 to 62% Brix measuring range
- **Model VBR62**
Handheld Refractometer with 28~62% Brix measuring range, 0.2% resolution
- **Model VHN2**
Handheld Refractometer with 58~90% Brix measuring range, 0.2% resolution
- **Model VSA1**
Handheld Refractometer with 0~100‰ / 1.000~1.070 Salinity measuring range and 1 ‰ / 0.005 resolution