

FH(P) series Fume Hood

Advantage

1. It is safer to use anti-corrosive water tap
2. Microprocessor control system, LED display
3. With memory function in case of power-failure
4. Made of porcelain white PP, resistant to strong acid, alkali and anti-corrosion.
5. Front window which is made of thick transparent toughened glass maximize light and visibility inside the fume hood, providing a bright and open working environment.



SPECIFICATIONS :

Model		FHI 1000(P)	FHI 1200(P)	FHI 1500(P)	FHI 1800(P)
External Size (WxDxH)		1040x750x2450 mm	1240x800x2450 mm	1540x800x2450 mm	1840x800x2450 mm
Internal Size (WxDxH)		820x520x872 mm	980x570x700 mm	1280 570 700 mm	1580x570x700 mm
Work Surface Height		820mm			
Max Opening		740mm			
Air Velocity		0.3~0.8m/s			
Noise		≤68dB			
Illuminating Lamp		LED lamp			
		12Wx1	30Wx2	30Wx2	36Wx2
Blower		Built-in PP centrifugal blower (2 blowers for FHI 1800(P) only) Speed adjustable			
Front Window		Resitant to acid and alkali, Manual, 5mm toughened glass, height adjustable.			
Power Supply		AC220V±10%, 50/60Hz; 110V±10%, 60Hz			
Consumption		330W	360W	360W	360W
Material	Main Body	Made of Porcelain white PP, thickness 8mm, resistant to strong acid, alkali and anti-corrosion			
	Work Table	Chemical resistant phenolic resin			
Standard Accessory		Illuminating Lamp, Water tap, Gas tap, Water Sink, Base cabinet			
		Waterproof socketx2, PP centrifugal blower, pipe strapx2 (4 Pcs for FH1800(P) only)			
		4 meters PVC duct (2 Pcs of 4 meters PVC dust for FHI 1800(P) only), Diameter : 250 mm			
Optional Accessory		PP work table, epoxy resin board or ceramic board.			
		Outer PVC centrifugal blower (only need to select when the duct more than 4 m)			
Gross Weight		160 kg	198 kg	225 kg	259 kg
Package Size (WxDxH)	Main Body	1190x890x1620 mm	1390x940x1620 mm	1690x940x162 mm	1990x940x1620 mm
	Base Cabinet	1190x890x1040 mm	1390x940x1040 mm	1690x940x1040mm	1990x940x1040 mm