



*Laboratory Equipment*

## **TDS Meter EC1385**

- Measuring range : 0.00~19.99EC, CF: 0.0~199CF, TDS: 10~19990ppm
- Accuracy :  $\pm 2\%$  F·S
- Resolution EC: 0.01EC, CF: 0.1CF, TDS: 10ppm
- Operating Temperature : 0~50°C
- Temperature Compensation : 0~50°C
- Power Supply 4×1.5V (AG13)
- Dimensions : 185×38×38mm



***INNOTECH Inc***

2648E Workman Ave # 413

West Covina, CA 91791

Tel. : 626-420759

Email : [innotech@innotechlab.com](mailto:innotech@innotechlab.com)

<http://www.innotechlab.com>