

Electric-Heated Thermostatic Water Bath



Features:

1. Internal material:
High quality stainless steel
2. Use microcomputer intellective temperature control instrument P.I.D, control temperature stablely, high precision, digital display, watch directly and clearly

Model	SY-1L2H	SY-2L4H	SY-2L6H	SY-2L8H
Temperature range	Chamber temperature +5 °C - 100 °C			
Temperature fluctuation	±0.5 °C			
Consumption power	500W	1000W	1500W	2000W
Power	110V-220V / 50-60Hz			
Chamber size (mm)	300 x 155 x 145	300 x 300 x 145	445 x 300 x 145	670 x 365 x 165
Physical dimension (mm)	Single line double holes	Double Lines four holes	Double Lines six holes	Double Lines eight holes



INNOTECH Inc.

2648E Workman Ave # 413

West Covina, CA 91791

Tel : 626-420756

Email : innotech@innotechlab.com