

# BiochemisBiochemistry Incubator



## Features :

1. Working area material: Stainless steel or galvanization. Round arc inner angle easy to clean
2. LED display. Touch screen set and adjust
3. Audio and visual alarm
4. Time range : 0-99 hours and 0-89 mins or no time limit
5. Double-layer door structure: open outer door, observe studio experimental results throught the inner door which is made from high intensity toughened glass, temperature and humidity unaffected
6. Adjustable shelves

Model	BJPX-Diego	BJPX-Hartford	BJPX-Dover	BJPX-Newark	BJPX-Warren
Capacity	60 L	80 L	150 L	250 L	300 L
Temperature range	0°C-60°C				
Temperature control	0.1°C				
Temperature fluctuation	±0.5°C				
Shelves	2	2	3	3	3
Time Range	0-99 hours and 0-59 mins or cintinuous				
Consumption power	250W	350W	450W	700W	800W
Power	110V-220V / 50-60Hz				
Cooling	R12 (normal) R134a (environment protection type)				
Chamber size (mm)	380 x 300 x 440	440 x 320 x 520	480 x 390 x 780	570 x 500 x 850	570 x 540 x 950
External dimension (mm)	510 x 525 x 970	570 x 545 x 1100	610 x 620 x 1435	700 x 735 x 1500	700 x 770 x 1600



**INNOTECH Inc.**

2648E Workman Ave # 413

West Covina, CA 91791

Tel : 626-420756

Email : [innotech@innotechlab.com](mailto:innotech@innotechlab.com)